



INDIA'S
1 BRAND
Experience the Reliable Technology
User-friendly Best Quality
"Confident"
Dental & Medical Equipments

MANUAL ICU COT
CME- ICUM-046

FEATURES :

- Height adjustable, manually operated, four section cots.
- Soft U Foam bed, 4 inch thick covered with Vinyl.
- 4 inch thick memory foam bed with vinyl cover is also provided, Which adjusts to body contour for free air circulation to prevent bed sore (optional).
- Head and footboards are easily detachable, made up of SS material single piece, can be cleaned by any cleaning solution.
- Bed Top: 4 section of MS perforated sheet.
- Mounted on 5 inches size 4 wheels with diagonal locking system.
- Provision for Saline stand on all 4 corners.
- Height, backrest, knee rest and trendelenberg operated thro four screw systems with thrust bearings individually manoeuvred by a single handle
- Backrest and Kneerest shall retract as they are individually raised
- Back rest operating systems shall have a load balancing mechanism.
- Castor housing wheels are made from high grade non floor-staining synthetic materials with integrated thread guards.
- ICU bed of minimum 130Kg weight bearing capacity
- SS Components wherever used of 304 grade



Technical Data

Variation up to $\pm 10\%$

• Length (Buffer to Buffer)	:	2000 mm
• Width (Buffer to Buffer)	:	900 mm
• Length (Top flatbed)	:	2100 mm
• Width (Top flatbed)	:	800 mm
• Height (w/o foam)	:	520 mm
• Backrest	:	0 - 40°
• Fowler's movement	:	0 - 30°



ISO 9001:2015
www.tuv.com
ID: 9105015614



CONFIDENT DENTAL EQUIPMENTS PVT. LTD.

Plot No. 17-H, Sector I, Phase II, Bidadi Industrial Area,
Bidadi - 562109, Bangalore (INDIA) GSTIN : 29AAACC9198C1ZF
Ph: +91-80 28023000 | E-mail: info@confidental.org | Web: www.confidental.org